

# BIOMEDICAL RESEARCH SERVICE CENTER

## UNIVERSITY at BUFFALO, STATE UNIVERSITY of NEW YORK

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 3435 Main Street, Buffalo, NY 14214, USA  
Tel/Fax: (716) 829-3106 Email: chunglee@buffalo.edu Web: www.bmrservice.com

### Aldehyde Dehydrogenase (ALDH) Assay Kit (Cat #: E-112 NAD<sup>+</sup> formula)

**COMPONENTS:** ALDH Assay Solution- 10 ml (for 200 assays); store in aliquots at -80°C after first thawing  
10x ALDH Substrate- 1 ml, store at -80°C  
10x Cell Lysis Solution- 25 ml, store at 4°C

**PRODUCT DESCRIPTION:** The assay solution uses acetaldehyde as a substrate and thus preferentially detects the activity of ALDH2. This assay is based on the reduction of INT to INT-formazan, which exhibits an absorption maximum at 492 nm ( $\epsilon = 18 \text{ mM}^{-1}\text{cm}^{-1}$ ) and allows for sensitive detection of ALDH activity in serum and tissue samples. Assay solution is stable for several years if stored and handled properly.

#### Preparation of cell/tissue extracts:

1. Wash  $\sim 10^6$  cells with ice-cold phosphate-buffered saline (PBS). Animal tissue should be washed with PBS thoroughly to remove blood cells. Freeze cell pellet and tissue at -80°C until use.
2. Add 50 – 100  $\mu\text{l}$  ice-cold 1x Cell Lysis Solution (diluted 10-fold with dH<sub>2</sub>O from 10x Cell Lysis Solution) to cell pellet. Extract cells by pipetting up and down (gently but thoroughly). Leave lysate on ice for 5 min with intermittent gentle agitation. Centrifuge lysate in a refrigerated microfuge for 3 min at maximum speed ( $\sim 13,000$  rpm) and harvest supernatant for ALDH assay. For tissue extraction, weigh  $\sim 50$  mg tissue and homogenize in 1 ml ice-cold 1x Cell Lysis Solution. Centrifuge homogenate at 4°C for 3 min at maximum speed and harvest supernatant. Store cell lysate and tissue homogenate at -80°C.
3. Perform protein assay to determine sample protein concentration. Equalize sample protein concentration by diluting with ice-cold 1x Cell Lysis Solution. A suggested protein concentration range is 0.5 – 2 mg/ml.

#### Reagent thawing:

Keep thawed ALDH Assay Solution and 10x ALDH Substrate on ice. Gently agitate solution prior to pipetting. It is important to minimize the time the reagents are thawed. Freeze solutions immediately after use.

#### Preparation of control solution and reaction solution:

Control solution is prepared by mixing 1 part of dH<sub>2</sub>O and 10 parts of ALDH Assay Solution, e.g. 50  $\mu\text{l}$  dH<sub>2</sub>O mixed with 500  $\mu\text{l}$  ALDH Assay Solution. Keep freshly prepared control solution on ice during assay.

Reaction solution is prepared by mixing 1 part of 10x ALDH substrate and 10 parts of ALDH Assay Solution, e.g. 50  $\mu\text{l}$  10x ALDH Substrate mixed with 500  $\mu\text{l}$  ALDH Assay Solution. Keep freshly prepared reaction solution on ice during assay. **Warning:** The ALDH substrate contains acetaldehyde and should be handled with caution in a fume hood.

Estimate the volume of control solution and reaction solution required for each assay. Each sample is treated with 50  $\mu\text{l}$  control solution and 50  $\mu\text{l}$  reaction solution in separate wells (see below).

#### Enzyme assay:

1. Add 10  $\mu\text{l}$  of each sample to a 96-well plate in duplicate: one set for control and another set for reaction.
2. After all samples have been pipetted to the plate in duplicate, add 50  $\mu\text{l}$  control solution to one set of wells and 50  $\mu\text{l}$  reaction solution to another set of wells. Gently agitate plate for 30 sec. Cover plate and incubate in a humidified 37°C incubator for 30 min or 60 min (for low activity). Do not use CO<sub>2</sub> incubator.
3. Stop assay by adding 50  $\mu\text{l}$  3% Acetic acid (not included in the kit) to each control solution well and reaction solution well followed by brief gentle agitation. Measure O.D.<sub>492 nm</sub> using a plate reader. Subtract control well reading from reaction well reading for each sample. Use the subtracted sample reading (**O.D.**) for enzyme activity calculation shown below.
4. Sample ALDH activity in IU/L unit =  $\mu\text{mol}/(\text{L}\cdot\text{min}) = \text{O.D.} \times 1000 \times 110 \mu\text{l} / (30 \text{ min} \times 0.6 \text{ cm} \times 18 \times 10 \mu\text{l}) = \text{O.D.} \times 33.96$ . If incubation for 60 min, ALDH activity in IU/L unit = **O.D. × 16.98**. Note that sample dilution may be desired to achieve assay linearity. Multiply the result by the dilution factor where applicable.

#### Additional information:

- A 3% Acetic acid solution needs to be prepared for reaction termination.
- The “0.6 cm” in the equation is the typical light path in a 96-well plate and may be custom adjusted as needed.
- The assay solution contains DMSO and iodinitrotetrazolium violet. The ALDH substrate contains acetaldehyde. Please refer to the product page of our website or contact us for MSDS information.