

# BIOMEDICAL RESEARCH SERVICE CENTER

## UNIVERSITY at BUFFALO, STATE UNIVERSITY of NEW YORK

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 3435 Main Street, Buffalo, NY 14214, USA  
Tel/Fax: (716) 829-3106 Email: chunglee@buffalo.edu Web: www.bmrservice.com

### Alkaline Phosphatase (ALP) Assay Kit (Cat #: E-102)

**COMPONENTS:** ALP Assay Solution- 40 ml (for 400 assays); **store in aliquots at -80°C after first thawing**  
10x Cell Lysis Solution- 25 ml, store at 4°C (**swirl bottle briefly prior to pipetting**)

**PRODUCT DESCRIPTION:** Alkaline phosphatase (ALP or ALK) possesses a hydrolytic activity removing phosphate groups from protein and non-protein molecules. The assay kit is based on conversion of p-nitrophenol phosphate to nitrophenol in an alkaline buffer. Nitrophenol exhibits an absorption maximum at 405 nm (molar extinction coefficient =  $18.75 \text{ mM}^{-1}\text{cm}^{-1}$ ), allowing for sensitive detection of ALP activity in crude cell/tissue extracts, serum/plasma and urine. Assay solution is stable for several years if stored and handled properly.

**Plasma:** Plasma samples may need to be diluted with ice-cold 1x Cell Lysis Solution (from 10-fold dilution of 10x Cell Lysis Solution with  $\text{dH}_2\text{O}$ ) to obtain assay linearity. Samples should be store at -80°C.

#### Preparation of cell/tissue extracts:

1. Wash  $10^5 - 10^6$  cells with ice-cold phosphate-buffered saline (PBS). Animal tissue should be washed with PBS thoroughly to remove blood cells. Freeze cell pellet and tissue at -80°C until use.
2. Add 50 – 100  $\mu\text{l}$  ice-cold 1x Cell Lysis Solution to cell pellet. Extract cells by pipetting up and down (gently but thoroughly). Leave lysate on ice for 5 min with intermittent gentle agitation. For tissue extraction, weigh ~25 mg tissue and homogenize in 0.5 ml ice-cold 1x Cell Lysis Solution. Briefly spin homogenate at ~2,000 rpm for 5 sec in a microfuge to deposit insoluble tissue debris. The lysate can be subjected to a high-speed centrifugation to separate the membrane and cytosolic protein fractions if necessary. Store lysate at -80°C.
3. Perform protein assay to determine sample protein concentration. Equalize sample protein concentration by diluting with ice-cold 1x Cell Lysis Solution. A suggested protein concentration range is 0.2 – 1 mg/ml.

**Enzyme assay for clear sample:** Thaw ALP Assay Solution and keep solution on ice prior to assay.

1. Add 10  $\mu\text{l}$  sample to each well of a 96-well plate. Add 10  $\mu\text{l}$   $\text{dH}_2\text{O}$  to a well serving as blank.
2. Reaction is initiated by adding 100  $\mu\text{l}$  ALP Assay Solution to each well. Mix contents by brief gentle agitation. Cover plate and incubate at 37°C for 30 min or 60 min.
3. Stop reaction by adding 20  $\mu\text{l}$  1N NaOH (not included in the kit) to each well followed by brief gentle agitation. Yellowish reaction color appears in wells during incubation.
4. Measure  $\text{O.D.}_{405 \text{ nm}}$  using a plate reader. Subtract blank well reading from each sample well reading to generate **O.D.**. For 30 min incubation, ALP enzyme activity in IU/L unit =  $\mu\text{mol}/(\text{L}\cdot\text{min}) = \text{O.D.} \times 1000 \times 130 \mu\text{l}/(30 \text{ min} \times 0.6 \text{ cm} \times 18.75 \times 10 \mu\text{l}) = \text{O.D.} \times 38.52$ . For 60 min incubation, ALP activity in IU/L unit = **O.D. x 19.26**.

#### Enzyme assay for hemolyzed sample:

1. Add 10  $\mu\text{l}$  sample to each well of a 96-well plate in duplicate.
2. Add 100  $\mu\text{l}$   $\text{dH}_2\text{O}$  to one set of sample wells. Add 100  $\mu\text{l}$  ALP Assay Solution to the other set of sample wells. Mix contents by brief agitation. Cover plate and incubate at 37°C for 30 min or 60 min.
3. Stop reaction by adding 20  $\mu\text{l}$  1N NaOH (not included in the kit) to each  $\text{dH}_2\text{O}$  well and assay well followed by brief agitation.
4. Measure  $\text{O.D.}_{405 \text{ nm}}$  using a plate reader. Subtract  $\text{dH}_2\text{O}$  well reading from assay well reading for each sample to generate **O.D.**. For 30 min incubation, ALP enzyme activity in IU/L unit =  $\mu\text{mol}/(\text{L}\cdot\text{min}) = \text{O.D.} \times 1000 \times 130 \mu\text{l}/(30 \text{ min} \times 0.6 \text{ cm} \times 18.75 \times 10 \mu\text{l}) = \text{O.D.} \times 38.52$ . For 60 min incubation, ALP activity in IU/L unit = **O.D. x 19.26**.

#### Additional information:

- Note that “0.6 cm” in the equation is the typical light path in a 96-well plate and may be custom adjusted as needed.
- A 1N NaOH solution needs to be prepared for reaction termination. Avoid skin contact with NaOH. Please refer to the product page of our website or contact us for MSDS information.