

BIOMEDICAL RESEARCH SERVICE CENTER

UNIVERSITY at BUFFALO, STATE UNIVERSITY of NEW YORK

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 3435 Main Street, Buffalo, NY 14214, USA

Tel/Fax: (716) 829-3106

Email: chunglee@buffalo.edu

Web: www.bmrservice.com

β -Hydroxybutyrate Dehydrogenase (BDH) Assay Kit (Cat #: E-139)

COMPONENTS: BDH Assay Solution- 5 ml (for 100 wells), store at -70°C (**shield solution from light during assay**)
10x BDH Substrate- 0.5 ml, store at -70°C
10x Cell Lysis Solution- 25 ml, store at 4°C (**contains 1% TX-100; swirl bottle briefly prior to dilution**)

PRODUCT DESCRIPTION: The BDH enzyme activity assay is based on the reduction of INT in a NADH-coupled reaction to formazan, which exhibits an absorption maximum at 492 nm ($\epsilon = 18 \text{ mM}^{-1}\text{cm}^{-1}$) and allows for sensitive measurement of BDH activity in tissue/cell extracts. Reagents are stable for several years if stored and handled properly.

Preparation of cell/tissue extracts:

1. Prepare 1x Cell Lysis Solution by diluting 10x Cell Lysis Solution with ice-cold dH₂O. Bring up at least $\sim 10^5$ washed cells in 50 – 100 μl ice-cold 1x Cell Lysis Solution by pipetting up and down gently. Leave lysate on ice for 5 min with agitation. If lysate is overly turbid, add more 1x Cell Lysis Solution and repeat pipetting. Tissue is homogenized in ice-cold 1x Cell Lysis Solution (10 – 20 mg tissue in 0.5 ml).
2. Centrifuge lysate in a cold microfuge at $\sim 14,000$ rpm for 5 min. Supernatant is harvested and stored at -80°C.
3. Use the BCA protein assay method to determine lysate protein concentration. A suggested sample protein concentration range is 1 – 2 mg/ml. Keep lysates on ice at all times during assay.

Reagent thawing:

Keep thawed BDH Assay Solution and 10x BDH Substrate on ice shielded from light. Do not over thaw. Gently agitate assay solution prior to first pipetting. It is important to minimize the time the reagents are thawed. Freeze solutions immediately after use.

Preparation of control solution and reaction solution:

Control solution is prepared by mixing 1 part of dH₂O and 9 parts of BDH Assay Solution, e.g. 50 μl dH₂O mixed with 450 μl BDH Assay Solution. Keep solution on ice.

Reaction solution is prepared by mixing 1 part of 10x BDH Substrate and 9 parts of BDH Assay Solution, e.g. 50 μl 10x BDH Substrate mixed with 450 μl BDH Assay Solution. Keep solution on ice and use immediately.

Enzyme assay:

1. Thaw lysates quickly and keep on ice (do not over thaw). Each sample is treated with 50 μl control solution and 50 μl reaction solution. Add 20 μl of each sample to a plain (uncoated) 96-well plate placed on ice in duplicate.
2. After all samples have been pipetted to the plate in duplicate, swiftly add 50 μl control solution to one set of wells and 50 μl reaction solution to the other set of wells. Mix contents by gentle agitation for 10 sec. Cover plate and incubate in a 37°C incubator for 30 – 60 min (do not use CO₂ incubator). Cherry red color should gradually appear in wells.
3. Measure O.D._{492 nm} using a plate reader at 30 min and at 60 min.
4. Subtract control well reading from reaction well reading for each sample for each time point. Use the subtracted reading ($\Delta\text{O.D.}$) for enzyme activity calculation. If incubation for 30 min, sample BDH activity in IU/L = $\mu\text{mol}/(\text{L}\cdot\text{min}) = \Delta\text{O.D.} \times 1000 \times 70 \mu\text{l} / (30 \text{ min} \times 0.5 \text{ cm} \times 18 \times 20 \mu\text{l}) = \Delta\text{O.D.} \times 12.96$. If incubation for 60 min, BDH activity = $\Delta\text{O.D.} \times 6.48$. Enzyme activity can be presented as units/ μg proteins.
Note: Incubation time may be increased to 120 min to increase $\Delta\text{O.D.}$ for samples exhibiting low BDH activity.

Additional information:

- The assay solution contains DMSO and iodinitrotetrazolium violet. Please refer to the product page of our website or contact us for MSDS information.