

BIOMEDICAL RESEARCH SERVICE, UNIVERSITY AT BUFFALO

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 955 Main St., Suite 4102, Buffalo, NY 14203
Tel: (716) 4586903 Email: chunglee@buffalo.edu Web: www.bmrservice.com

11 β -Hydroxysteroid dehydrogenase type 2 (11 β -HSD2) Assay Kit (Cat #: E-152)

COMPONENTS: 11 β -HSD2 Assay Solution- 5 ml (for 100 wells), store at -80°C (**shield solution from light during assay**)
10x 11 β -HSD2 Control- 0.5 ml, store at 4°C
10x 11 β -HSD2 Reaction- 0.5 ml, store at 4°C
10x Cell Lysis Solution- 25 ml, store at room temperature (**agitate vial to redissolve contents upon arrival**)

PRODUCT DESCRIPTION: The 11 β -HSD2 enzyme activity assay is based on the oxidation of cortisol (corticosterone), which is coupled to NADH-dependent reduction of INT to formazan. The formazan product exhibits an absorption maximum at 492 nm (extinction coefficient = 18 mM⁻¹cm⁻¹), allowing for measurement of 11 β -HSD2 activity present in cell/tissue extracts and serum. Kit components are stable for at least one year if stored and handled properly.

Preparation of cell/tissue extracts:

1. Prepare 1x Cell Lysis Solution by diluting 10x Cell Lysis Solution with dH₂O. Bring up ~10⁵ washed cells in 50 – 100 μ l ice-cold 1x Cell Lysis Solution by pipetting up and down gently. Leave lysate on ice for 5 min with agitation. If lysate is overly turbid, add more 1x Cell Lysis Solution and repeat pipetting. Tissue is homogenized in ice-cold 1x Cell Lysis Solution (~5 mg tissue in 0.5 ml).
2. Centrifuge lysate in a cold microfuge at ~14,000 rpm for 5 min. Supernatant is harvested and stored at -80°C.
3. Use the BCA protein assay method to determine lysate protein concentration. A suggested sample protein concentration range is 0.5 – 4 mg/ml. *Normalization of lysate protein concentration is recommended.*

Reagent thawing:

11 β -HSD2 Assay Solution: Keep thawed assay solution on ice and shielded from light. *Do not over thaw.* Freeze assay solution immediately after use.

10x 11 β -HSD2 Control and Reaction: Keep both solutions at room temperature prior to assay. Vortex solution to dissolve precipitates if necessary.

Preparation of control solution and reaction solution:

Control solution is prepared by mixing 1 part of 10x 11 β -HSD2 Control and 10 parts of 11 β -HSD2 Assay Solution, e.g. 50 μ l 11 β -HSD2 Control mixed with 500 μ l 11 β -HSD2 Assay Solution.

Reaction solution is prepared by mixing 1 part of 10x 11 β -HSD2 Reaction and 10 parts of 11 β -HSD2 Assay Solution, e.g. 50 μ l 11 β -HSD2 Reaction mixed with 500 μ l 11 β -HSD2 Assay Solution.

Enzyme assay:

1. Keep thawed samples on ice. *Do not over thaw.* Add 20 μ l of each cell/tissue lysate to duplicate wells of a plain (uncoated) 96-well plate.
2. Add 50 μ l control solution to one set of wells and 50 μ l reaction solution to the other set of wells. Mix contents by gentle agitation for 10 sec. Cover plate and incubate in a 37°C humidified incubator for 60 min. *Do not use CO₂ incubator.*
3. Measure O.D._{492 nm} using a plate reader. Subtract control well reading from reaction well reading for each sample. Use the subtracted reading (Δ O.D.) for enzyme activity calculation.
4. Sample 11 β -HSD2 enzyme activity in IU/L = μ mol/(L·min) = Δ O.D. $\times$ 1000 $\times$ 70 μ l / (60 min $\times$ 0.4 cm $\times$ 18 $\times$ 20 μ l) = Δ O.D. $\times$ 8.1. Enzyme activity may be normalized by GAPDH enzyme activity (#E-101) or protein concentration (#A-117).

Additional information:

- The assay solution contains dimethyl sulfoxide (DMSO) and iodonitrotetrazolium (INT). 10x 11 β -HSD2 Reaction contains cortisol. Please refer to the product webpage or contact us for material SDS information.
- 10x Cell Lysis Solution precipitates at cold temperature. Redissolve contents by agitation in warm water upon arrival.