

BIOMEDICAL RESEARCH SERVICE, UNIVERSITY AT BUFFALO

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 955 Main St., Suite 4102, Buffalo, NY 14203
Tel: (716) 8293106 Email: chunglee@buffalo.edu Web: www.bmrservice.com

Acid α -Glucosidase (GAA) Assay kit (Cat #: E-148)

COMPONENTS: GAA Assay Solution- 10 ml (for 100 wells), store at -20°C
GAA Control Solution- 10 ml, store at -20°C
10x Cell Lysis Solution- 25 ml, store at room temperature (**agitate vial to redissolve contents**)

PRODUCT DESCRIPTION: The lysosomal acid α -glucosidase (GAA or acid maltase) assay is based on the hydrolysis of 4-nitrophenyl- α -D-glucopyranoside. The enzyme is capable of degrading glucose- α -1,4-glucose and glucose- α -1,6-glucose glycosidic bonds. Cleavage of nitrophenyl from the substrate increases absorbance at 405 nm (extinction coefficient = 18 mM⁻¹cm⁻¹), allowing for measurement of GAA enzyme activity present in tissue/cell lysates and serum. Kit components are stable for at least 1 year if stored and handled properly.

Preparation of cell/tissue extracts:

1. Prepare 1x Cell Lysis Solution by diluting 10x Cell Lysis Solution with ice-cold dH₂O. Bring up ~10⁵ washed cells in 100 – 200 μ l ice-cold 1x Cell Lysis Solution by pipetting up and down gently. Leave lysate on ice for 5 min with agitation. If lysate is overly turbid, add more 1x Cell Lysis Solution and repeat pipetting. Tissue is homogenized in ice-cold 1x Cell Lysis Solution (~10 mg tissue in 0.5 ml).
2. Centrifuge lysate in a cold microfuge at ~14,000 rpm for 5 min. Supernatant is harvested and stored at -70°C.
3. Use the BCA protein assay method to determine lysate protein concentration. A suggested sample protein concentration range is 0.2 – 2 mg/ml.

Enzyme assay for cell/tissue lysates:

1. Thaw Assay Solution and keep solution on ice. *Do not over thaw*. Slight precipitates may appear after thawing. Agitate container to dissolve contents if necessary. Add 20 μ l of each cell/tissue sample to a plain (uncoated) 96-well plate. Add 20 μ l 1x Cell Lysis Solution to an empty well serving as blank.
2. Reaction is initiated by adding 100 μ l GAA Assay Solution to each well. Mix contents by brief gentle agitation. Cover plate and incubate at 37°C for 60 min (do not use CO₂ incubator).
3. Measure O.D._{405 nm} using a plate reader. Subtract blank well reading from each sample well reading to generate Δ O.D..

Enzyme assay for serum:

1. Thaw Assay Solution and Control Solution and keep solutions on ice. *Do not over thaw*. Slight precipitates may appear after thawing. Vortex container to dissolve contents if necessary. Add 5 μ l of each serum sample to duplicate wells of a plain (uncoated) 96-well plate (for control set and reaction set).
2. Add 100 μ l GAA Control Solution to one set of sample wells (control). Add 100 μ l GAA Assay Solution to the other set of sample wells (reaction). Mix contents by brief gentle agitation. Cover plate and incubate at 37°C for 60 min (do not use CO₂ incubator).
3. Measure O.D._{405 nm} using a plate reader. Subtract control well reading from reaction well reading to generate Δ O.D. for each sample.

Enzyme activity calculation:

For 20 μ l sample, GAA enzyme activity in IU/L unit = μ mol/(L•min) = Δ O.D. x 1000 x 120 μ l / (60 min x 0.6 cm x 18 x 20 μ l) = Δ O.D. x 9.26. For 5 μ l sample, GAA enzyme activity in IU/L unit = μ mol/(L•min) = Δ O.D. x 1000 x 105 μ l / (60 min x 0.6 cm x 18 x 5 μ l) = Δ O.D. x 32.41. Enzyme activity may be presented as units/ μ g proteins or normalized by GAPDH enzyme activity (see cat#E-101).

Additional information:

- Please refer to the product page of our website or contact us for SDS information.
- 10x Cell Lysis Solution precipitates upon freezing. Re-dissolve contents by agitation in warm water upon arrival.