

BIOMEDICAL RESEARCH SERVICE, UNIVERSITY AT BUFFALO

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 955 Main St., Suite 4102, Buffalo, NY 14203
Tel: (716) 4586903 Email: chunglee@buffalo.edu Web: www.bmrservice.com

Cathepsin B (CTSB) Assay kit (Cat #: E-146)

COMPONENTS: 10x CTSB Buffer- 0.6 ml, store at room temperature
100x CTSB Substrate- 50 μ l, store at -20°C (for 100 wells)
100x DTT: 0.1 ml, store at -20°C
10x Cell Lysis Solution- 25 ml, store at room temperature (**agitate vial to redissolve contents**)

PRODUCT DESCRIPTION: Cathepsin B (CTSB) is a lysosomal cysteine protease of the papain family. The CTSB activity assay kit is based on proteolytic hydrolysis of the chromogenic peptide substrate z-Arg-Arg-p-nitroanilide. Cleavage of nitroanilide from the substrate increases absorbance at 405 nm (extinction coefficient= 9.96 mM⁻¹cm⁻¹), allowing for measurement of CTSB activity present in tissue/cell lysates and biological fluids such as serum/plasma. Kit components are stable for several years if stored and handled properly.

Preparation of cell/tissue extracts:

1. Prepare 1x Cell Lysis Solution by diluting 10x Cell Lysis Solution with ice-cold dH₂O. Bring up ~10⁵ washed cells in 100 – 200 μ l ice-cold 1x Cell Lysis Solution by pipetting up and down gently. Leave lysate on ice for 5 min with agitation. If lysate is overly turbid, add more 1x Cell Lysis Solution and repeat pipetting. Tissue is homogenized in ice-cold 1x Cell Lysis Solution (~10 mg tissue in 0.5 ml).
2. Centrifuge lysate in a cold microfuge at ~14,000 rpm for 5 min. Supernatant is harvested and stored at -80°C.
3. Use the BCA protein assay method to determine lysate protein concentration. A suggested sample protein concentration range is 0.2 – 2 mg/ml.

Preparation of working solution:

Each well of a 96-well plate requires 50 μ l of freshly prepared working solution. For preparation of 1 ml of working solution, for example, add 100 μ l 10x CTSB Buffer to 900 μ l ddH₂O followed by addition of 10 μ l 100x DTT and 10 μ l 100x CTSB Substrate. Mix contents well and use the working solution immediately.

Enzyme assay:

1. Calculate the amount of working solution required for each assay, and prepare the solution as described above.
2. Place a plain (uncoated) 96-well plate on one sheet of Kimwipe positioned on ice. Add 10 μ l of each sample to the plate.
3. Add 50 μ l working solution to each sample well. Mix contents by brief gentle agitation for 10 sec. Blot dry bottom of plate. IMMEDIATELY measure O.D._{405 nm} using a plate reader. Data are recorded as O.D._{0 min}.
4. Cover plate and incubate at 37°C for 30 and 60 min. *Do not use CO₂ incubator.* Measure O.D._{405 nm} after incubation. Data are recorded as O.D._{30 min} and O.D._{60 min}.
5. If incubation for 30 min, CTSB activity in IU/L = (O.D._{30 min} – O.D._{0 min}) x 1000 x 60 μ l / (30 min x 0.3 cm x 9.96 x 10 μ l) = (O.D._{30 min} – O.D._{0 min}) x **66.93**. If incubation for 60 min, CTSB activity in IU/L = (O.D._{60 min} – O.D._{0 min}) x **33.47**

Enzyme activity can be normalized by protein concentration and presented as units/ μ g proteins. Alternatively, enzyme activity may be normalized by GAPDH enzyme activity (see cat#E-101).

Additional information:

- The 100x CTSB Substrate solution contains dimethyl sulfoxide (DMSO). Please refer to the product webpage or contact us for SDS information.
- 10x Cell Lysis Solution precipitates at cold temperature. Re-dissolve contents by agitation in warm water upon arrival.