

BIOMEDICAL RESEARCH SERVICE, UNIVERSITY AT BUFFALO

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 955 Main St., Suite 4102, Buffalo, NY 14203
Tel: (716) 4586903 Email: chunglee@buffalo.edu Web: www.bmrservice.com

Glutamate Dehydrogenase (GDH or GLDH) Assay Kit (Cat #: E-123)

COMPONENTS: GDH Assay Solution- 5 ml (for 100 wells); store at -80°C (**shield solution from light during assay**)
10x GDH Substrate- 0.5 ml, store at -20°C or -80°C
10x Cell Lysis Solution- 25 ml, store at room temperature (**agitate vial to redissolve contents**)

PRODUCT DESCRIPTION: The GDH enzyme activity assay is based on the reduction of the tetrazolium salt INT in a NADH-coupled enzymatic reaction to formazan, which exhibits an absorption maximum at 492 nm ($\epsilon = 18 \text{ mM}^{-1}\text{cm}^{-1}$) and allows for measurement of GDH activity in tissue/cell extracts and serum/plasma. Kit components are stable for several years if stored and handled properly.

Preparation of cell/tissue extracts:

1. Prepare 1x Cell Lysis Solution by diluting 10x Cell Lysis Solution with dH₂O. Bring up $\sim 10^5$ washed cells in 50–100 μl ice-cold 1x Cell Lysis Solution by pipetting up and down gently. Leave lysate on ice for 5 min with agitation. If lysate is overly turbid, add more 1x Cell Lysis Solution and repeat pipetting. Tissue is homogenized in ice-cold 1x Cell Lysis Solution ($\sim 5 \text{ mg}$ tissue in 0.5 ml).
2. Centrifuge lysate in a cold microfuge at $\sim 14,000 \text{ rpm}$ for 5 min. Supernatant is harvested and stored at -80°C.
3. Use the BCA protein assay method to determine lysate protein concentration. A suggested sample protein concentration range is 0.1 – 3 mg/ml. *Normalization of lysate protein concentration is recommended.*

Reagent thawing:

Keep thawed GDH Assay Solution and 10x GDH Substrate on ice. *Do not over thaw.* Freeze solutions immediately after use.

Preparation of control solution and reaction solution:

Control solution is prepared by mixing 1 part of dH₂O and 10 parts of GDH Assay Solution, e.g. 50 μl dH₂O + 500 μl GDH Assay Solution. Keep solution on ice.

Reaction solution is prepared by mixing 1 part of 10x GDH Substrate and 10 parts of GDH Assay Solution, e.g. 50 μl 10x GDH Substrate + 500 μl GDH Assay Solution. Keep solution on ice.

Enzyme assay:

1. Add 10 μl of each sample to duplicate wells of a plain (uncoated) 96-well plate.
2. Add 50 μl control solution to one set of wells (control wells) and 50 μl reaction solution to the other set of wells (reaction wells). Mix contents by gentle agitation for 10 sec. Cover plate and incubate in a 37°C incubator for 30 min or 60 min. *Do not use CO₂ incubator.* Cherry red color should gradually appear in sample wells.
3. Measure O.D._{492 nm} using a plate reader.
4. Subtract control well reading from reaction well reading for each sample. Use the subtracted reading ($\Delta\text{O.D.}$) for enzyme activity calculation.
5. If incubation for 30 min, GDH enzyme activity in IU/L = $\mu\text{mol}/(\text{L}\cdot\text{min}) = \Delta\text{O.D.} \times 1000 \times 60 \mu\text{l} / (30 \text{ min} \times 0.4 \text{ cm} \times 18 \times 10 \mu\text{l}) = \Delta\text{O.D.} \times 27.78$. If incubation for 60 min, GDH activity = $\Delta\text{O.D.} \times 13.89$.
6. Enzyme activity may be presented as units/ μg proteins or normalized by GAPDH enzyme activity (see cat#E-101). Increase protein concentration or incubation time if detected enzyme activity is too low.

Additional information:

- The assay solution contains dimethyl sulfoxide (DMSO) and iodonitrotetrazolium (INT). Please refer to the product webpage or contact us for material SDS information.
- 10x Cell Lysis Solution precipitates at cold temperature. Redissolve contents by agitation in warm water upon arrival.