

BIOMEDICAL RESEARCH SERVICE, UNIVERSITY AT BUFFALO

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 955 Main St., Suite 4102, Buffalo, NY 14203
Tel: (716) 4586903 Email: chunglee@buffalo.edu Web: www.bmrservice.com

Glutaminase (GLS) Assay Kit (Cat #: E-133)

COMPONENTS: TA Assay Solution- 5 ml (for 100 wells); store at -80°C (**shield solution from light during assay**)
GLS Buffer- 5 ml, store at -80°C
Glutamine- 50 mg, store at 4°C
10x Sample Buffer- 25 ml, store at 4°C

PRODUCT DESCRIPTION: The glutaminase (GLS) enzyme activity assay is based on sequential GLS-mediated hydrolytic reaction and glutamate dehydrogenase-mediated redox reaction, the latter of which is coupled to reduction of INT to formazan (molar extinction coefficient $\epsilon = 18 \text{ mM}^{-1}\text{cm}^{-1}$ at 492 nm), allowing for measurement of GLS activity present in crude cell/tissue lysates. Kit components are stable for several years if handled and stored properly.

Preparation of cell/tissue extracts:

1. PBS-washed cell/tissue samples should be frozen at -80°C prior to homogenization. At least $\sim 10^5$ cells or 2 mg of tissue should be used for lysate preparation. Freeze thawing, mechanical grinding, or sonication facilitates protein extraction.
2. Prepare 1x Sample Buffer by diluting 10x Sample Buffer with dH₂O. Bring up freeze-thawed cells in 50 – 100 μl ice-cold 1x Sample Buffer by pipetting up and down. Adherent cells after deep freezing can be directly harvested from plate by scraping. Leave lysate on ice for 5 min with agitation. Freeze-thawed tissue is mechanically homogenized in ice-cold 1x Sample Buffer ($\sim 5 \text{ mg}$ tissue in 500 μl).
3. Centrifuge lysate in a cold microfuge at $\sim 14,000 \text{ rpm}$ for 5 min. Supernatant is harvested and stored at -80°C.
4. Perform protein assay to determine lysate (supernatant) protein concentration. A suggested sample protein concentration range is 0.2 – 2 mg/ml. *Normalization of sample protein concentration is recommended.*

Preparation of glutamine solution: Assay of each sample requires 40 μl glutamine solution. Prepare a fresh 5 mg/ml glutamine solution in GLS Buffer, e.g., 2.5 mg glutamine in 0.5 ml GLS Buffer. Vigorously vortex solution to dissolve glutamine. Discard unused glutamine solution.

TA Assay Solution: Keep thawed TA Assay Solution on ice and shielded from light. *Do not over thaw*. Freeze solution immediately after use.

Enzyme assay:

1. Use a plain (uncoated) 96-well plate for the assay. Pipette 10 μl of each sample to duplicate wells (for control and reaction).
2. Add 40 μl GLS Buffer to each control well. Add 40 μl Glutamine solution to each reaction well. Gently agitate plate for 10 sec to mix. Cover plate and incubate in a humidified 37°C incubator for 30 min. *Do not use CO₂ incubator.*
3. Remove plate from incubator. Add 50 μl TA Assay Solution to each control well and reaction well. Gently agitate plate to mix. Cover plate and incubate in a humidified 37°C incubator for another 30 min.
4. Remove plate from incubator. Measure O.D.492 nm using a plate reader.
5. Subtract control well reading from reaction well reading for each sample. Use the subtracted reading ($\Delta\text{O.D.}$) for enzyme activity calculation as shown below.
6. $\text{GLS enzyme activity in IU/L} = \mu\text{mol}/(\text{L} \cdot \text{min}) = (\Delta\text{O.D.} \times 1000 \times 100 \mu\text{l}) \div (30 \text{ min} \times 0.4 \text{ cm} \times 18 \times 10 \mu\text{l}) = \Delta\text{O.D.} \times 46.3$. Enzyme activity may be presented as units/ μg proteins or normalized by GAPDH enzyme activity (Cat#E-101).

Additional information:

The kit contains dimethyl sulfoxide (DMSO) and iodonitrotetrazolium (INT). Please refer to the product webpage or contact us for material SDS information.