

## **BIOMEDICAL RESEARCH SERVICE, UNIVERSITY AT BUFFALO**

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 955 Main St., Suite 4102, Buffalo, NY 14203  
Tel: (716) 4586903 Email: chunglee@buffalo.edu Web: www.bmrservice.com

### **Succinate Dehydrogenase (SDH) Assay Kit (Cat #: E-161)**

**COMPONENTS:** SDH Assay Solution- 5 ml (for 100 wells); store at -80°C (**shield solution from light during assay**)  
20x SDH Substrate- 0.25 ml, store at -80°C  
10x Lysis Buffer- 25 ml, store at 4°C (dilute 10-fold with ice-cold dH<sub>2</sub>O)  
10x Extraction Buffer- 20 ml, store at 4°C (dilute 10-fold with ice-cold dH<sub>2</sub>O)

**PRODUCT DESCRIPTION:** The SDH enzyme activity assay is based on reduction of the tetrazolium salt INT in a FADH<sub>2</sub>-coupled reaction to formazan, which exhibits an absorption maximum at 492 nm ( $\epsilon = 18 \text{ mM}^{-1}\text{cm}^{-1}$ ), allowing for measurement of SDH activity present in cell/tissue extracts. Kit components are stable for several years if stored and handled properly.

#### **Preparation of cell/tissue lysates:**

1. **Cells:** Bring up  $\sim 10^6$  freeze-thawed cells in  $\sim 0.1$  ml ice-cold 1x Lysis Buffer by pipetting up and down. Centrifuge lysate in a cold microfuge at  $\sim 14,000$  rpm for 5 min. Collect supernatant for assays such as Hexokinase, GAPDH, and LDH. Resuspend pellet (membrane fraction) in  $\sim 50$   $\mu\text{l}$  ice-cold 1x Extraction Buffer by pipetting up and down thoroughly, making sure the pellet is homogenized. Leave homogenate on ice for  $\sim 10$  min with agitation. Repeat manual homogenization. Centrifuge homogenate at  $\sim 5000$  rpm for 1 min. Harvest supernatant for SDH assay. Recommended sample protein concentration is 0.2 - 2 mg/ml.
2. **Tissue:** Mechanically homogenize freeze-thawed tissue in ice-cold 1x Lysis Buffer ( $\sim 5$  mg tissue in 0.5 ml). Centrifuge lysate in a cold microfuge at  $\sim 14,000$  rpm for 5 min. Collect supernatant for assays such as Hexokinase, GAPDH, and LDH. Resuspend pellet (membrane fraction) in  $\sim 0.25$  ml ice-cold 1x Extraction Buffer by pipetting up and down thoroughly, making sure the pellet is homogenized. Leave homogenate on ice for  $\sim 10$  min with agitation. Repeat manual homogenization. Centrifuge homogenate at  $\sim 5000$  rpm for 1 min. Harvest supernatant for SDH assay. Recommended sample protein concentration is 0.2 - 2 mg/ml.

#### **Preparation of control solution and reaction solution**

Thaw SDH Assay Solution and 20x SDH Substrate and keep solutions on ice shielded from light. *Do not over thaw*. Prepare fresh control and assay solutions as follows:

Control solution is prepared by mixing 1 part of dH<sub>2</sub>O and 20 parts of SDH Assay Solution, e.g. 25  $\mu\text{l}$  dH<sub>2</sub>O + 500  $\mu\text{l}$  SDH Assay Solution.

Reaction solution is prepared by mixing 1 part of 20x SDH Substrate and 20 parts of SDH Assay Solution, e.g. 25  $\mu\text{l}$  20x SDH Substrate + 500  $\mu\text{l}$  SDH Assay Solution.

#### **Enzyme assay**

1. Thaw protein samples and keep on ice. *Do not over thaw*. Pipette each sample up and down several times prior to sample addition. Add 20  $\mu\text{l}$  of each sample to duplicate wells of a plain (uncoated) 96-well plate (for control set and reaction set).
2. Add 50  $\mu\text{l}$  control solution to one set of sample wells and 50  $\mu\text{l}$  reaction solution to the other set of sample wells. Gently agitate plate for 10 sec to mix contents. Cover plate and incubate in a 37°C incubator for 60 min or 90 min. *Do not use CO<sub>2</sub> incubator*.
3. Measure O.D.<sub>492 nm</sub> using a plate reader. Subtract control well reading from reaction well reading for each sample. Use the subtracted reading ( $\Delta\text{O.D.}$ ) for enzyme activity calculation.
4. If incubation for 60 min, SDH activity in IU/L =  $\mu\text{mol}/(\text{L}\cdot\text{min}) = \Delta\text{O.D.} \times 1000 \times 70 \mu\text{l}/(60 \text{ min} \times 0.4 \text{ cm} \times 18 \times 20 \mu\text{l}) = \Delta\text{O.D.} \times 8.1$  (or  $\Delta\text{O.D.} \times 5.4$  for 90 min). SDH activity can be normalized by protein concentration or GAPDH activity (#E-101).

#### **Additional information:**

- The assay solution contains dimethyl sulfoxide (DMSO) and iodonitrotetrazolium (INT). Please refer to the product webpage or contact us for material SDS information.