

BIOMEDICAL RESEARCH SERVICE, UNIVERSITY AT BUFFALO

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 955 Main St., Suite 4102, Buffalo, NY 14203
Tel: (716) 4586903 Email: chunglee@buffalo.edu Web: www.bmrservice.com

Threonine Dehydrogenase (TDH) Assay Kit (Cat #: E-132)

COMPONENTS: TDH Assay Solution- 5 ml (for 100 wells), store at -80°C (**shield solution from light during assay**)
10x TDH Substrate- 0.5 ml, store at -80°C
10x Sample Buffer- 25 ml, store at 4°C

PRODUCT DESCRIPTION: The TDH enzyme activity assay is based on reduction of the tetrazolium salt INT in NADH-coupled enzymatic reaction to formazan, which exhibits an absorption max at 492 nm (molar extinction coefficient = $18 \text{ mM}^{-1}\text{cm}^{-1}$) and allows for measurement of TDH activity in crude cell/tissue extracts. Kit components are stable for several years if stored and handled properly.

Preparation of cell/tissue extracts:

1. PBS-washed cell/tissue samples should be frozen at -80°C prior to homogenization. At least $\sim 10^5$ cells or 2 mg of tissue should be used for lysate preparation. Freeze thawing, mechanical grinding, or sonication facilitates protein extraction.
2. Prepare 1x Sample Buffer by diluting 10x Sample Buffer with dH₂O. Bring up freeze-thawed cells in 50–100 μl ice-cold 1x Sample Buffer by pipetting up and down. Adherent cells after deep freezing can be directly harvested from plate by scraping. Leave lysate on ice for 5 min with agitation. Freeze-thawed tissue is mechanically homogenized in ice-cold 1x Sample Buffer ($\sim 5 \text{ mg}$ tissue in 500 μl).
3. Centrifuge lysate in a cold microfuge at $\sim 14,000 \text{ rpm}$ for 5 min. Supernatant is harvested and stored at -80°C.
4. Perform protein assay to determine lysate (supernatant) protein concentration. A suggested sample protein concentration range is 0.5 – 3 mg/ml. *Normalization of sample protein concentration is recommended.*

Reagent thawing:

Keep thawed TDH Assay Solution and 10x TDH Substrate on ice. *Do not over thaw*. Keep TDH Assay Solution shielded from light. Freeze solutions immediately after use.

Preparation of control solution and reaction solution:

Control solution is prepared by mixing 1 part of dH₂O and 10 parts of TDH Assay Solution, e.g. 50 μl dH₂O mixed with 500 μl TDH Assay Solution.

Reaction solution is prepared by mixing 1 part of 10x TDH Substrate and 10 parts of TDH Assay Solution, e.g. 50 μl 10x TDH Substrate mixed with 500 μl TDH Assay Solution.

Enzyme assay:

1. Use a plain (uncoated) 96-well for the assay. Add 20 μl of each sample to duplicate wells (for control set and reaction set).
2. Add 50 μl control solution to one set of sample wells (control) and 50 μl reaction solution to the other set of sample wells (reaction). Mix contents by gentle agitation for 10 sec. Cover plate and incubate in a 37°C incubator for 60 min or 120 min. *Do not use CO₂ incubator*. Cherry red color should gradually appear in wells.
3. Measure O.D._{492 nm} using a plate reader. Subtract control well reading from reaction well reading for each sample. Use the subtracted reading ($\Delta\text{O.D.}$) for enzyme activity calculation.
4. If incubation for 60 min, TDH enzyme activity in $\text{IU/L} = \mu\text{mol}/(\text{L}\cdot\text{min}) = \Delta\text{O.D.} \times 1000 \times 70 \mu\text{l}/(60 \text{ min} \times 0.4 \text{ cm} \times 18 \times 20 \mu\text{l}) = \Delta\text{O.D.} \times 8.1$ (or $\Delta\text{O.D.} \times 4.05$ for 120 min). *Increase sample protein concentration if low enzyme activity is observed*. Enzyme activity may be presented as units/ μg proteins or normalized by GAPDH enzyme activity (see cat#E-101).

Additional information:

The assay solution contains dimethyl sulfoxide (DMSO) and iodonitrotetrazolium (INT). Please refer to the product webpage or contact us for material SDS information.