

BIOMEDICAL RESEARCH SERVICE, UNIVERSITY AT BUFFALO

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 955 Main St., Suite 4102, Buffalo, NY 14203
Tel: (716) 8293106 Email: chunglee@buffalo.edu Web: www.bmrservice.com

ATP Assay Kit (Cat #: A-107)

COMPONENTS: ATP Assay Solution: 2 x 10 ml, store at -80°C (for 200 assays); **shield solution from light during assay**
1 mM (1000 μ M) ATP Standard: 0.1 ml, store at -80°C

PRODUCT DESCRIPTION: The enzyme luciferase catalyzes, in the presence of ATP, the oxidation of luciferin with concomitant emission of yellow-green light (560 nm), which can be conveniently measured by a luminometer (Anal. Biochem. 80:496,1977). This provides the basis for a sensitive non-radioactive ATP assay system capable of detecting 0.1 μ M or 1 picomole of ATP. The ATP Assay Solution if stored in aliquots at -80°C and shielded from light is stable for many years. Repeated freeze-thaw cycles of ATP Assay Solution should be avoided.

PROTOCOLS

ATP Assay Solution: Each ATP measurement requires 0.1 ml ATP Assay Solution. Store solution in aliquots after the first thawing. Thawed ATP Assay Solution should be kept on ice shielded from light. Do not over-thaw.

ATP Standard: Dilute the 1 mM (1000 μ M) ATP Standard with dH₂O to 0.5 – 10 μ M (see graph below).

Sample deproteinization:

The TCA method (*Methods in Enzymology* 133:14-22,1986) is recommended for cell and tissue samples. Please follow the protocol at <https://bmrservice.com/deproteinization-by-tca>. Store TCA- and ether-treated samples at -20°C. TCA-treated samples may be diluted 3-fold with dH₂O prior to assay to achieve optimal result.

Attention: TCA is highly corrosive. Ether is volatile and flammable. Exercise caution handling the reagents.

The boiling water method (*Anal Biochem* 306:323-327, 2002) may be suitable for some cell samples. Pellet $\sim 10^6$ saline-washed cells in a microtube, and remove saline completely. Add 100 μ l ice-cold dH₂O and disrupt the cell pellet by vigorous vortexing. Immediately heat the lysate in a boiling water bath for 10 min, following which the boiled lysate is clarified by centrifugation at $\sim 14,000$ rpm) for 5 min. Store supernatant at -80°C.

ASSAY

1. Thaw enough ATP Assay Solution (0.1 ml required for each measurement) and keep solution on ice shielded from light.
2. Add 10 μ l of each ATP standard and sample to a clear (uncoated) 96-well plate. Add 10 μ l dH₂O to a well for background Relative Light Unit (RLU) reading.
3. Add 100 μ l ATP Assay Solution to each well. Immediately measure light emission using a luminometer.
4. The background (H₂O) RLU reading should be subtracted from all standard and sample RLU readings. Use subtracted RLU to calculate sample ATP concentration.
5. Generate a plot of ATP standard RLU vs. ATP concentrations with trendline equation. A new plot must be generated for each assay.
6. Calculate sample ATP concentration using the derived equation (x = sample ATP concentration in μ M; y = sample RLU).

