

BIOMEDICAL RESEARCH SERVICE, UNIVERSITY AT BUFFALO

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 955 Main St., Suite 4102, Buffalo, NY 14203
Tel: (716) 4586903 Email: chunglee@buffalo.edu Web: www.bmrservice.com

Glutathione (GSH) Assay Kit (Cat #: A-126)

COMPONENTS:
GSH Assay Solution: 50 ml, store at 4°C
Glutathione (GSH): 100 mg, store at 4°C
200x CDNB: 0.1 ml, store at -80°C
50x GST: 0.4 ml, store at -80°C (100 assays)

PRODUCT DESCRIPTION: Glutathione (GSH) is a tripeptide, γ -L-glutamyl-L-cysteinyl-glycine, involved in diverse cellular function. The GSH assay kit is based on a spectrophotometric assay (PNAS 71:3879-3882,1974), which employs 1-chloro-2,4-dinitrobenzene (CDNB) as a substrate for GST-mediated GSH conjugation. The conjugation causes an increase in the absorbance at 340 nm. The assay requires a UV-transparent 96-well plate or 0.2-ml quartz cuvette for measuring changes in O.D._{340 nm}. Assay reagents are stable for several years if stored and handled properly.

PROTOCOLS

Sample preparation by TCA extraction:

Cell/tissue samples are deproteinized by TCA and extracted by diethyl ether prior to assay. The extraction protocol is available at <https://bmrservice.com/deproteination-by-tca>. *Use caution handling TCA and ether.*

Glutathione (GSH) standard:

First prepare a 4 mM GSH (5 mg GSH dissolved in 4 ml dH₂O). Dilute the solution 10-fold with dH₂O to obtain 0.4 mM (400 μ M) GSH. Perform 1:1 serial dilution with dH₂O to obtain 200, 100, and 50 μ M GSH standards. GSH standards should be prepared fresh for each assay.

Reagent thawing:

Keep thawed 200x CDNB at room temperature. Briefly vortex the CDNB vial to homogenize solution. Keep thawed 50x GST solution on ice (*do not over thaw GST solution*). Freeze the reagents immediately after use.

Preparation of working solution:

Each GSH standard and sample requires 0.2 ml working solution. Estimate the volume of working solution required for each experiment and prepare fresh. An example for 10 reactions: 2 ml GSH Assay Solution + 10 μ l 200x CDNB. Vortex solution immediately to prevent CDNB precipitation. Proceed to add 40 μ l 50x GST to the solution followed by brief mixing. Keep working solution on ice shielded from light and use immediately.

ASSAY

1. Add 20 μ l dH₂O (blank), 20 μ l of each GSH standard and 20 μ l of each sample to a UV-transparent 96-well plate. Alternatively, use a 0.2-ml quartz cuvette for the UV measurement.
2. Add 0.2 ml freshly prepared working solution to each well or tube. Mix contents by gentle agitation for 5 sec. Incubate at 37°C shielded from light for 20 min.
3. Measure O.D._{340 nm} using a plate reader or spectrophotometer. Presence of GSH will cause increased O.D._{340 nm} reading in a concentration-dependent manner.
4. Subtract blank reading from each GSH standard and sample reading to generate respective Δ O.D.. Generate a graph of GSH concentrations (x axis) vs. Δ O.D. (y axis) and a trendline equation.
5. Apply sample Δ O.D. to the equation to obtain sample GSH concentration. A new GSH standard plot must be generated for each assay.

Additional information:

- 200x CDNB contains dimethyl sulfoxide (DMSO) solvent and 1-chloro-2,4-dinitrobenzene. Please refer to the product page of our website or contact us for material SDS information.

