

BIOMEDICAL RESEARCH SERVICE, UNIVERSITY AT BUFFALO

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 955 Main St., Suite 4102, Buffalo, NY 14203

Tel: (716) 4586903

Email: chunglee@buffalo.edu

Web: www.bmrservice.com

Short-Chain Fatty Acid Oxidation (FAO) Assay Kit (Cat #: E-141S)

COMPONENTS: FAO Assay Solution- 5 ml (for 100 wells), store at -80°C (**shield solution from light during assay**)
25x FAO Substrate (S)- 0.2 ml, store at -80°C (**avoid multiple freeze-thaw cycles**)
10x Sample Buffer- 25 ml, store at 4°C

PRODUCT DESCRIPTION: The FAO assay is based on oxidation of butyryl-CoA coupled to FADH₂/NADH-dependent reduction of idonitrotetrazolium to formazan, which exhibits an absorption maximum at 492 nm ($\epsilon = 18 \text{ mM}^{-1}\text{cm}^{-1}$). The kit is for measurement of FAO activity in cell/tissue lysates. Reagents are stable for at least one year if stored and handled properly.

Preparation of cell/tissue extracts:

1. PBS-washed cell/tissue samples should be frozen at -80°C prior to homogenization. At least $\sim 10^5$ cells or 2 mg of tissue should be used for lysate preparation. Freeze thawing, mechanical grinding, or sonication facilitates protein extraction.
2. Prepare 1x Sample Buffer by diluting 10x Sample Buffer with dH₂O. Bring up freeze-thawed cells in 50 – 100 μl ice-cold 1x Sample Buffer by pipetting up and down. Adherent cells after deep freezing can be directly harvested from plate by scraping. Leave lysate on ice for 5 min with agitation. Freeze-thawed tissue is mechanically homogenized in ice-cold 1x Sample Buffer (~ 2 mg tissue in 100 μl).
3. Centrifuge lysate in a cold microfuge at $\sim 14,000$ rpm for 5 min. Supernatant is harvested and stored at -80°C.
4. Perform protein assay to determine lysate (supernatant) protein concentration. A suggested sample protein concentration range is 0.5 – 4 mg/ml. *Normalization of sample protein concentration is recommended.*

Reagent thawing:

Keep thawed FAO Assay Solution and 25x FAO Substrate (S) on ice and shielded from light. *Do not over thaw*. Freeze solutions immediately after use.

Preparation of control solution and reaction solution:

Control solution is prepared by mixing 1 part of dH₂O and 25 parts of FAO Assay Solution, e.g. 20 μl dH₂O + 500 μl FAO Assay Solution.

Reaction solution is prepared by mixing 1 part of 25x FAO Substrate (S) and 25 parts of FAO Assay Solution, e.g. 20 μl 25x FAO Substrate (S) + 500 μl FAO Assay Solution.

Enzyme assay:

1. Thaw lysates quickly and keep on ice. *Do not over thaw*. Add 20 μl of each sample to duplicate wells of a plain (uncoated) 96-well plate. For drug discovery application, mix 1 μl of an inhibitor with sample in wells first.
2. Add 50 μl control solution to one set of sample wells and 50 μl reaction solution to the other set of sample wells. Mix contents by gentle agitation for 10 sec. Cover plate and incubate in a 37°C humidified incubator for 60-120 min. *Do not use CO₂ incubator*. Cherry red color should gradually appear in wells.
3. Remove plate from incubator. Eliminate air bubbles from wells if necessary. Measure absorbance (O.D.) at 492 nm using a plate reader.
4. Subtract control well reading from reaction well reading for each sample. Use the subtracted reading ($\Delta\text{O.D.}$) for enzyme activity calculation. If incubation for 60 min, sample FAO activity in $\text{IU/L} = \mu\text{mol}/(\text{L}\cdot\text{min}) = \Delta\text{O.D.} \times 1000 \times 70 \mu\text{l} / (60 \text{ min} \times 0.4 \text{ cm} \times 18 \times 20 \mu\text{l}) = \Delta\text{O.D.} \times 8.10$. If incubation for 120 min, FAO activity = $\Delta\text{O.D.} \times 4.05$. Enzyme activity may be normalized by GAPDH enzyme activity (see cat#E-101) or protein concentration.

Additional information:

- The assay solution contains dimethyl sulfoxide (DMSO) and idonitrotetrazolium (INT). Please refer to the product webpage or contact us for material SDS information.